



Recombinant Arabidopsis thaliana Cysteine-rich repeat secretory protein 57 (CRRSP57)

Product Code	CSB-MP866928DOA
Abbreviation	CRRSP57
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FHD5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	ISESDSDEHM ATFCNDSSGN FTRNTTYNTN LNTLLSTLSN QSSFANYYNL TTGLGSDTVH GMFLCIGDVN RTTCNACVKN ATIEIAKNCT NHREAIYYF SCMVRYSDFK FLSTLETKPN TYWSSDDPIP KSYDKFGQRL SDKMGVEIIR SSLLSSSFTP YYLMDTTTFD NLYDLESVVQ CSPHLDPKNC TTCLKLALQE LTQCCGDQLW AFIFTPKCLV SFDTSNSSSL PPLPPPSRSG S
Source	Mammalian cell
Target Names	CRRSP57
Protein Names	Recommended name: Cysteine-rich repeat secretory protein 57
Expression Region	21-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.