



Recombinant Arabidopsis thaliana Putative FBD-associated F-box protein At5g22720 (At5g22720)

Product Code	CSB-YP866997DOA
Abbreviation	At5g22720
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FNI9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEETKVKRTC LERTVCSSKA KEDLISQLPD SLITQILFYI QTKKAVTTSV LSKRWRSLWL STPGLVLISN DFTDYNFVSV FVDKFLGFSR EQKLCLHKLK LSIRKGENDQ DCVTRWIDFV ATPKPKHLDV EIGPTRCECF EVIPLSLYSC ESLLYLRLNH VCLGKFESVS LPCLKTMSLE QNIYANEADL ESLISTCPVL EDLSFVSGAY DKVNVLRVQS QTLTSLNIEG CVEYLDLTKS EVLIDATRLK YLNLEADQYE SKTIRNSGSL TKVNULLGYFH MKNNDDDDEV DLQKRDMVHN FFTSISGVSD MKISSQAFAL FIMNTMPFSP KFCNLSCLEV EIFLPSLETL PTFLESFPNL KSLILGLRYW TPKKELRLSF VPRCLLSSLE FVEIKGCSRS NMERVKYVGE PIETKLARYF VENSTILKKL VLPLRSSTPE EIYSVDFWNF LEFPRRSSIC QITYVAGL
Source	Yeast
Target Names	At5g22720
Protein Names	Recommended name: Putative FBD-associated F-box protein At5g22720
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.