



Recombinant Arabidopsis thaliana Protein NIM1-INTERACTING 1 (NIMIN-1)

Product Code	CSB-EP867001DOA
Abbreviation	NIMIN-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FNZ5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MYPKQFSLYN YSLETMSKDE NVESKETIRV DKRVREDEEE EEEKKIDTFF KLIKHYQEAR KRRREELAEN SGVRRKSNV GERSGIVVPA FQPEDFSQCR TGLKPPLMFV SDHKEENTKV EQEEDQTEER NEDKALDLNL AL
Source	E.coli
Target Names	NIMIN-1
Protein Names	Recommended name: Protein NIM1-INTERACTING 1 Short name= Protein NIMIN-1
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.