



Recombinant Arabidopsis thaliana Transcription factor GTE6 (GTE6)

Product Code	CSB-EP867012DOA-B
Abbreviation	GTE6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FT54
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MADSVPGHVA GGGLQGFSDV AECIKQRVDE VLQWVDSLEH KLKEVEEFYS SIGVSNVSGSI GKDTEKGRHV VGIRKIQQEA ARREAVAAKR MQDLMRQFGT IFRQITQHKC AWPFMHPVNV EGLGLHDYFE VIDKPMDFST IKNQMEAKDG TGYKHVMQIY ADMRLVFENA MNYNEETSDV YSMAKKLLEK FEEKWAHFLP KVQEEEEKIRE EEEKQAAKEA LLAKEASHIK TTRELGNEIC HANDELEKLM RKVVERCRKI TIEEKRNIGL ALLKLSPDDL QKVLGIVAQA NPSFQPRAEE VSIEMDILDE PTLWRLKFFV KDALDNAMKK KKEEETKTRE LSGAQKKEVS KKRNATTKLA ERKTKRSRI
Source	E.coli
Target Names	GTE6
Protein Names	Recommended name: Transcription factor GTE6 Alternative name(s): Bromodomain-containing protein GTE6 Protein GENERAL TRANSCRIPTION FACTOR GROUP E6 Protein GLOBAL TRANSCRIPTION FACTOR GROUP E6
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.