



# Recombinant Arabidopsis thaliana Transcription factor GTE6 (GTE6)

<b>Product Code</b>	CSB-YP867012DOA
<b>Abbreviation</b>	GTE6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FT54
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADSVPGHVA GGGLQGFSVD AECIKQRVDE VLQWVDSLEH KLKEVEEFYS SIGVSNVSGSI GKDTEKGRHV VGIRKIQQEA ARREAVAAKR MQDLMRQFGT IFRQITQHKC AWPFMHPVNV EGLGLHDYFE VIDKPMDFST IKNQMEAKDG TGYKHVMQIY ADMRLVFENA MNYNEETSDV YSMAKKLLEK FEEKWAHFLP KVQEEEEKIRE EEEKQAAKEA LLAKEASHIK TTRELGNEIC HANDELEKLM RKVVERCRKI TIEEKRNIGL ALLKLSPDDL QKVLGIVAQA NPSFQPRAEE VSIEMDILDE PTLWRLKFFV KDALDNAMKK KKEEETKTRE LSGAQKKEVS KKRNATTKLA ERKTKRSRI
<b>Source</b>	Yeast
<b>Target Names</b>	GTE6
<b>Protein Names</b>	Recommended name: Transcription factor GTE6 Alternative name(s): Bromodomain-containing protein GTE6 Protein GENERAL TRANSCRIPTION FACTOR GROUP E6 Protein GLOBAL TRANSCRIPTION FACTOR GROUP E6
<b>Expression Region</b>	1-369
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.