



Recombinant Protein male abnormal 21 (mab-21)

Product Code	CSB-EP867097CXX-B
Abbreviation	mab-21
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9GQ38
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	>85% (SDS-PAGE)
Sequence	MLGHNQNVVY QVNTFYNEKV QQRKARVTKN VHRIAKVVQE ILKEVEAQEP RFINTLTETS TGRFDGIVVH SPSEYEAVLY LNQMGVFNFV DDGTIQGCAV LKLSDGRKRS MSLWVEFITA SGYLSARKIR HRFQNIVAQV LQTPQFSEYC KLLQDNTDVR VRVDDKYTVQ ITCAFRNGI WPRSASHWPL AGLPWPNTAL ANQTKAEGFD LTSRETAITQ HNNPNKQAST MEADAWAMKM HGAENMLLTG GRQKTLSEILK CLRDAHMDFP GTPVTNYILK TLVLYECEKH CSEYEWEDTN IGDRLVGVLL QLVSCLQCRR CAHYFLPSLD LLRAKPTHTI EHSAKLTWHL VRKLMIDPNA LQTL
Source	E.coli
Target Names	mab-21
Protein Names	Recommended name: Protein male abnormal 21
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>