



Recombinant Human Uncharacterized protein C21orf56 (C21orf56)

Product Code	CSB-MP867117HU
Abbreviation	C21orf56
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0A9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEGGELMSR LLENADLKK QVRLLENQM LRLLSQSCQ EGGGHDLLPP RAHAYPEAGS PGSGVPDFGR FTSVADTPSQ LQTSSLEDLL CSHAPLSSD DTSPGCAAPS QAPFKAFLSP PEPHSHRGTD RKLSPLLSPL QDSLVDKTL EPREMVRPKK VCFSESSLPT GDRTRRSYYL NEIQSFAGAE KDARVVGEIA FQLDRRILAY VFPGVTRLYG FTVANIPEKI EQTSTKSLDG SVDERKLREL TQRYLALSAR LEKLGYSRDV HPAFSEFLIN TYGILKQRPD LRANPLHSSP AALRKLVIDV VPPKFLGDSL LLLNCLCELS KEDGKPLFAW
Source	Mammalian cell
Target Names	SPATC1L
Protein Names	Recommended name: Uncharacterized protein C21orf56
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.