



Recombinant Human Sharpin (SHARPIN)

Product Code	CSB-EP867119HU-B
Abbreviation	SHARPIN
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9H0F6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAPPAGGAAA AASDLGSAAV LLAVHAAVRP LGAGPDAAEQ LRRLQLSADP ERPGRFRLEL LGAGPGAVNL EWPLESVSYT IRGPTQHELQ PPPGGPGTLS LHFLNPQEAQ RWAVLVRGAT VEGQNGSKSN SPPALGPEAC PVSLPSPPEA STLKGPPPEA DLPRSPGNLT EREELAGSLA RAIAGGDEKG AAQVAAVLAQ HRVALSVQLQ EACFPPIR LQVTLEDAAS AASAASSAHV ALQVHPHCTV AALQEQVFSE LGFPPAVQRW VIGRCLCVPE RSLASYGVRQ DGDPAFLYLL SAPREAPATG PSPQHPQKMD GELGRLFPPS LGLPPGPQPA ASSLP SPLQP SWSCPSCTFI NAPDRPGCEM CSTQRPCTWD PLAAAST</p>
Source	E.coli
Target Names	SHARPIN
Protein Names	Recommended name: Sharpin Alternative name(s): Shank-associated RH domain-interacting protein Shank-interacting protein-like 1 Short name= hSIPL1
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>