



Recombinant Human Pseudouridylate synthase 7 homolog-like protein (PUS7L)

Product Code	CSB-YP867121HU
Abbreviation	PUS7L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0K6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEDTDYRIR FSSLCFFNDH VGFHGTIKSS PSDFIVIEID EQGQLVNKTI DEPIFKISEI QLEPNNFPKK PKLDLQNL SL EDGRNQEVHT LIKYTDGDQN HQSGSEKEDT IVDGTSKCEE KADVLSSFLD EKTHELLNNF ACDVREKWLS KTELIGLPPE FSIGRILDKN QRASLHSAIR QKFPFLVTVG KNSEIVVKPN LEYKELCHLV S EEEAFDFFK YLDAKKENSK FTFKPD TNKD HRKAVHHFVN KKFGNLVETK SFSKMNC SAG NPNVVVTVRF REKAHKRGKR PLSECQEGKV IYTAFTLRKE NLEMFEAIGF LAIKLGVIPS DFSYAGLKDK KAITYQAMVV RKVTPERLKN IEKEIEKKRM NVFNIRSVDD SLRLGQLKGN HFDIVIRNLK KQINDSANLR ERIMEAIENV KKKGFVNYYG PQRFGKGRKV HTDQIGLALL KNEMMKAIKL FLTPEDLDDP VNRACKYFLQ TEDAKGTL SL MPEFKVRERA LLEALHRFGM TEEGCIQAWF SLP HSMRIFY VHAYTSKIWN EAVSYRLETY GARVVQGD LV CLDEDID DEN FPNSKIHLVT EEEGSANMYA IHQVVL PVLG YNIQYPKNKV GQWYHDILSR DGLQTCRFKV PTLKLNIPGC YRQILKHPCN LSYQLMEDHD IDVKTKGSHI DETALSLLIS FDL DASCYAT VCLKEIMKHD V
Source	Yeast
Target Names	PUS7L
Protein Names	Recommended name: Pseudouridylate synthase 7 homolog-like protein EC=5.4.99.-
Expression Region	1-701
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.