



Recombinant Human Protein FAM113A (FAM113A)

Product Code	CSB-MP867128HU
Abbreviation	FAM113A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H1Q7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVFCLSSEEP RRPLRSDMVH FQASEVQQLL HNKFVVILGD SIQRAVYKDL VLLLQKDSLL TAAQLKAKGE LSFEQDQLVA GGQLGELHNG TQYREVRQFC SGSGHHLVRF YFLTRVYSEY LEDVLEELTY GPAPDLVIIN SCLWDLSRYG RCSMESYREN LERVFVRMDQ VLPDSCLLVW NMAMPLGERI TGGFLLPELQ PLAGSLRRDV VEGNFYSATL AGDHCDFVLD LHFHFRHAVQ HRHRDGVHWD QHAHRHLSHL LLTHVADAWG VELPKRGYPP DPWIEDWAEM NHPFQGSHRQ TPDFGEHLAL LPPPPSSLPP PMPFPYPLPQ PSPPPLFPPL PQDTPFFPGQ PFPPEFFNY NPVEDFSMPP HLGCGPGVNF VPGPLPPPPI GPNPHGQHWG PVVHRGMPRY VPNSPYHVRR MGGPCRQRLR HSERLIHTYK LDRRPPAHSG TWPG
Source	Mammalian cell
Target Names	PCED1A
Protein Names	Recommended name: Protein FAM113A Alternative name(s): Sarcoma antigen NY-SAR-23
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.