



# Recombinant Human Epsin-3 (EPN3)

<b>Product Code</b>	CSB-EP867130HU
<b>Abbreviation</b>	EPN3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H201
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTTSALRRQV KNIVHNYSEA EIKVREATSN DPWGPPSSLM SEIADLTFNT VAFTEVMGML WRRLNDSGKN WRHVYKALT LDYLLKTGSE RVAHQCRENL YTIQTLKDFQ YIDRDGKDQG VNVREKVKQV MALLKDEERL RQERTHALKT KERMALEGIG IGSGQLGFSR RYGEDYSRSR GSPSSYNSSS SSPRYTSDLE QARPQTSGEE ELQLQLALAM SREEAEKPVP PASHRDEDLQ LQLALRLSRQ EHEKEVRSWQ GDGSPMANGA GAVVHHQRDR EPEREERKEE EKLKTSQSSI LDLADIFVPA LAPPSTHCSA DPWDIPGFRP NTEASGSSWG PSADPWSPPI SGTVLSRSQP WDLTPMLSSS EPWGRTPVLP AGPPTTDPWA LNSPHHKLPS TGADPWGASL ETSDTPGGAS TFDPFKAPPE STETKEGLEQ ALPSGKPSSP VELDLFGDPS PSSKQNGTKE PDALDLGILG EALTQPSKEA RACRTPESFL GPSASSLVNL DSLVKAPQVA KTRNPFLTGL SAPSPTNPFGL AGEPGRPTLN QMRTGSPALG LAGGPVGAPL GSMTYSASLP LPLSSVPAGL TLPASVSVFP QAGAFAPQPL LPTPSSAGPR PPPPQTGTNP FL
<b>Source</b>	E.coli
<b>Target Names</b>	EPN3
<b>Protein Names</b>	Recommended name: Epsin-3 Alternative name(s): EPS-15-interacting protein 3
<b>Expression Region</b>	1-632
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.