



Recombinant Human BolA-like protein 2 (BOLA2)

Product Code	CSB-MP867143HU
Abbreviation	BOLA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H3K6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELSAEYLRE KLQRDLEAEH VEVEDTTLNR CSCSFRVLVV SAKFEGKPLL QRHRLVNAEL AEELPHIHAF EQKTLTPDQW ARERQK
Source	Mammalian cell
Target Names	BOLA2
Protein Names	Recommended name: BolA-like protein 2
Expression Region	1-86
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is located within a region of a segmental duplication on chromosome 16p. The product of this gene belongs to the eukaryotic subfamily of the BolA-like proteins. This gene encodes the BolA-like protein 2. The BolA-like proteins are widely conserved from prokaryotes to eukaryotes, and these proteins seem to be involved in cell proliferation or cell-cycle regulation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.