



Recombinant Human Pancreatic progenitor cell differentiation and proliferation factor (PPDPF)

Product Code	CSB-BP867146HU
Abbreviation	PPDPF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H3Y8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAIPSSGSL VATHDYRRR LGSTSSNSSC SSTECPGEAI PHPPGLPKAD PGHWWASFFF GKSTLPFMAT VLESAEHSEP PQASSSMTAC GLARDAPRKQ PGGQSSTASA GPPS
Source	Baculovirus
Target Names	PPDPF
Protein Names	Recommended name: Pancreatic progenitor cell differentiation and proliferation factor Alternative name(s): Exocrine differentiation and proliferation factor
Expression Region	1-114
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.