



# Recombinant Human GTP-binding protein 5 (GTPBP5)

<b>Product Code</b>	CSB-EP867150HU-B
<b>Abbreviation</b>	GTPBP5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H4K7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAPARCF SAR LRTVFQGVGH WALSTWAGLK PSRLLPQRAS PRLLSVGRAD LAKHQELPGK KLLSEKCLKR YFVDYRRVLV CGGNGGAGAS CFHSEPRKEF GGPDGGDGGN GGHVILRVDQ QVKSLSSVLS RYQGFSGEDG GSKNCFGRSG AVLYIRVPVG TLVKEGGRVV ADLSCVGD EY I AALGGAGGK GNRFFLANNN RAPVTCTPGQ PGQQRVLHLE LKTVAHAGMV GFPNAGKSSL LRAISNARPA VASYPFTTLK PHVGIVHYEG HLQIAVADIP GIIRGAHQNR GLGSAFLRHI ERCRFLLFVV DLSQPEPWTQ VDDLKYELEM YEKGLSARPH AIVANKIDL P EAQANLSQLR DHLGQEVIVL SALTGENLEQ LLLHLKVL YD AYAEAELGQG RQPLRW
<b>Source</b>	E.coli
<b>Target Names</b>	MTG2
<b>Protein Names</b>	Recommended name: GTP-binding protein 5 Alternative name(s): Protein obg homolog 1 Short name= ObgH1
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.