



Recombinant Human Ribosome production factor 2 homolog (RPF2)

Product Code	CSB-MP867164HU
Abbreviation	RPF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H7B2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDTLDRVVKP KTKRAKRFLE KREPKLNENI KNAMLIKGGN ANATVTKVLK DVYALKKPYG VLYKKNITR PFEDQTSLEF FSKKSDCSLF MFGSHNKKRP NNLVIGRMYD YHVLDMIELG IENFVSLKDI KNSKCPEGTK PMLIFAGDDF DVTEDYRRLK SLLIDFFRGP TVSNIRLAGL EYVLHFTALN GKIYFRSYKL LKKSGCRTP RIELEEMGPS LDLVLRRTL ASDDLYKLSM KMPKALKPKK KKNISHDTFG TTYGRIHMQK QDLSKLQTRK MKGLKKRPAE RITEDHEKKS KRIKKN
Source	Mammalian cell
Target Names	RPF2
Protein Names	Recommended name: Ribosome production factor 2 homolog Alternative name(s): Brix domain-containing protein 1 Ribosome biogenesis protein RPF2 homolog
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.