



Recombinant Human UPF0364 protein C6orf211 (C6orf211)

Product Code	CSB-MP867178HU
Abbreviation	C6orf211
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H993
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AVVPASLSG QDVGSFAYLT IKDRIPQILT KVIDTLHRHK SEFFEKHGEE GVEAEKKAIS LLSKLRNELQ TDKPFIPLVE KFDVTDIWNQ YLEYQQSLLN ESDGKSRWFY SPWLLVECYM YRRIHEAIQ SPPIDYFDVF KESKEQNFYG SQESIILCT HLQQLIRTIE DLDENQLKDE FFKLLQISLW GNKCDLSLSG GESSQNTNV LNSLEDLKP ILLNDMEHLW SLLSNCKKTR EKASATRVYI VLDNSGFELV TDLILADFL SSELATEVHF YGKTIPWFVS DTTIHDFNWL IEQVKHSNHK WMSKCGADWE EYIKMGKWVY HNHIFWTLPH EYCAMPQVAP DLYAELQKAH LILFKGDLNY RKL TGDRKWE FSVPFHQALN GFHPAPLCTI RTLKAEIQVG LQPGQGEQLL ASEPSWWTTG KYGIFQYDGP L
Source	Mammalian cell
Target Names	ARMT1
Protein Names	Recommended name: UPF0364 protein C6orf211
Expression Region	2-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.