



# Recombinant Human Putative E3 ubiquitin-protein ligase UNKL (UNKL)

<b>Product Code</b>	CSB-MP867181HU
<b>Abbreviation</b>	UNKL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H9P5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPSVSKAAAA ALSGSPPQTE KPTHYRYLKE FRTEQCPLFS QHKCAQHRPF TCFHWHLNQ RRRRPLRRRD GTFNYS PDVY CSKYNEATGV CPDGDECPYL HRTTGDTERK YHLRYKTGT CIHETDARGH CVKNGLHCAF AHGPLDLRPP VCDVRELQAQ EALQNGQLGG GEGVPDLQPG VLASQAMIEK ILSEDPRWQD ANFVLGSYKT EQCPKPPRLC RQGYACPHYH NSRDRRRNPR RFQYRSTPCP SVKHGDEWGE PSRCDGGDGC QYCHSRTEQQ FHPESTKCND MRQTGYCPRG PFCAFAHVEK SLGMVNEWGC HDLHLTSPSS TGSGQPGNAK RRDSPAEGGP RGSEQDSKQN HLAVFAAVHP PAPSVSSSSVA SSLASSAGSG SSSPTALPAP PARALPLGPA SSTVEAVLGS ALDLHLSNVN IASLEKDL EE QDGHDLGAAG PRSLAGSAPV AIPGSLPRAP SLHSPSSAST SPLGSLSQPL PGPVGSSAMT PPQQPPPLRS EPGTLGSAAS SYSPLGLNGV PGSIWDFVSG SFSPSPSPIL SAGPPSSSSA SPNGAELARV RRQLDEAKRK IRQWEESWQQ VKQVCDAWQR EAQEAKERAR VADSDRQLAL QKKEEVEAQV IFQLRAKQCV ACRERAHGAV LRPCQHHILC EPCAATAPEC PYCKGQLQW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UNKL
<b>Protein Names</b>	Recommended name: Putative E3 ubiquitin-protein ligase UNKL EC= 6.3.2.- Alternative name(s): RING finger protein unkempt-like Zinc finger CCCH domain-containing protein 5-like
<b>Expression Region</b>	1-680
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.