



Recombinant Human CaiB/baiF CoA-transferase family protein C7orf10 (C7orf10)

Product Code	CSB-BP867183HU
Abbreviation	C7orf10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAC7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> QS DMNNIKPLEG VKILDLTRVL AGPFATMNLG DLGAEVIKVE RPGAGDDTRT WGPPFVGTES TYYLSVNRNK KSIAVNIKDP KGVKIIKELA AVCDVFVENY VPGKLSAMGL GYEDIDEIAP HIIYCSITGY GQTGPISQRA GYDAVASAVS GLMHITGPEN GDPVRPGVAM TDLATGLYAY GAIMAGLIQK YKTGKGLFID CNLLSSQVAC LSHIAANYLI GQKEAKRWGT AHGSIVPYQA FKTKDGYIVV GAGNNQQFAT VCKILDPEL IDNSKYKTNH LRVHNRKELI KILSERFEEE LTSKWLYLFE GSGVPYGPIN NMKNVFAEPQ VLHNGLVMEM EHPTVGKISV PGPAVRYSKF KMSEARPPPL LGQHTTHILK EVLRYDDRAI GELLSAGVVD QHETH </p>
Source	Baculovirus
Target Names	SUGCT
Protein Names	Recommended name: CaiB/baiF CoA-transferase family protein C7orf10 EC= 2.-.-.- Alternative name(s): Dermal papilla-derived protein 13
Expression Region	39-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.