



# Recombinant Human CaiB/baiF CoA-transferase family protein C7orf10 (C7orf10)

<b>Product Code</b>	CSB-MP867183HU
<b>Abbreviation</b>	C7orf10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9HAC7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           QS DMNNIKPLEG VKILDLTRVL AGPFATMNLG DLGAEVIKVE            RPGAGDDTRT WGPPFVGTES TYYLSVNRNK KSIAVNIKDP KGVKIIKELA            AVCDVFVENY VPGKLSAMGL GYEDIDEIAP HIIYCSITGY GQTGPISQRA            GYDAVASAVS GLMHITGPEN GDPVRPGVAM TDLATGLYAY GAIMAGLIQK            YKTGKGLFID CNLLSSQVAC LSHIAANYLI GQKEAKRWGT AHGSIVPYQA            FKTKDGYIVV GAGNNQQFAT VCKILDPEL IDNSKYKTNH LRVHNRKELI            KILSERFEEE LTSKWLYLFE GSGVPYGPIN NMKNVFAEPQ VLHNGLVMEM            EHPTVGKISV PGPAVRYSKF KMSEARPPPL LGQHTTHILK EVLRYDDRAI            GELLSAGVVD QHETH         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SUGCT
<b>Protein Names</b>	Recommended name: CaiB/baiF CoA-transferase family protein C7orf10 EC= 2.-.-.- Alternative name(s): Dermal papilla-derived protein 13
<b>Expression Region</b>	39-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.