



Recombinant Human TLR adapter interacting with SLC15A4 on the lysosome (TASL)

Product Code	CSB-MP867184HU
Abbreviation	CXorf21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAI6
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLSEGYLSGL EYWNDIHWSC ASYNEQVAGE KEEETNSVAT LSYSSVDETQ VRSLYVSCKS SGKFISSVHS RESQHRSRQR VTVLQTNPNP VFESPNLAAV EICRDASRET YLVPSSCKSI CKNYNDLQIA GGQVMAINS SV TDFPSESSF EYGPLLKSS E IPLPMEDSIS TQPSPFPQKP IQRYSSYWR I TSIKEKSSLQ MQNPISNAVL NEYLEQKVVE LYKQYIMDTV FHDSSPTQIL ASELIMTSVD QISLQVSREK NLETSKARDI VFSRLLQLMS TEITEISTPS LHSQYSNVN P
Source	Mammalian cell
Target Names	CXorf21
Protein Names	Recommended name: Uncharacterized protein CXorf21
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.