



Recombinant Human RNA methyltransferase-like protein 1 (RNMTL1)

Product Code	CSB-MP867196HU
Abbreviation	RNMTL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HC36
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAALVRPARF VVRPLLQVVQ AWDLDARRWV RALRRSPVKV VFPSGEVVEQ KRAPGKQPRK APSEASAQEQ REKQPLEESA SRAPSTWEES GLRYDKAYPG DRRLSSVMTI VKSRPFREKQ GKILLEGRRL ISDALKAGAV PKMFFFSRLE YLKELPVDKL KGVSLIKVKF EDIKDWSDLV TPQGIMGIFA KPDHVKMTYP KTQLQHSLPL LLICDNLDP GNLGTILRSA AGAGCSKVLK TKGCVDWEP KVLRAGMGAH FRMPIINNLE WETVPNYLPP DTRVYVADNC GLYAQAEMSN KASDHGWVCD QRVMKFHKYE EEEDVETGAS QDWLPHVEVQ SYSDSWTEAP AAVVIGGETY GVSLESLQLA ESTGGKRLLI PVVPGVDSL SAMAASILLF EGKRQLRGRA EDLSRDRSYH
Source	Mammalian cell
Target Names	MRM3
Protein Names	Recommended name: RNA methyltransferase-like protein 1 EC= 2.1.1.-
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.