



Recombinant Human Zinc finger and BTB domain-containing protein 20 (ZBTB20)

Product Code	CSB-MP867198HU
Abbreviation	ZBTB20
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HC78
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLERKKPKTA ENQKASEENE ITQPGGSSAK PGLPCLNFEA VLSPDPALIH STHSLTNSHA HTGSSDCDIS CKGMTERIHS INLHNFSNSV LETLNEQRNR GHFCDVTVRI HGSMLRAHRC VLAAGSPFFQ DKLLLGYSI EIPSVVSVQS VQKLIDFMYG GVLRVVSQSEA LQILTAASIL QIKTVIDECT RIVSQNVGDV FPGIQDSGQD TPRGTPESGT SGQSSDTEG YLQSHQPQHSV DRIYSALYAC SMQNGSGERS FYSGAVVSHH ETALGLPRDH HMEDPSWITR IHERSQQMER YLSTTPETTH CRKQRPVRI QTLVGNIIHK QEMEDDYDYY GQQRVQILER NESEECTEDT DQAEGTESEP KGESFDSGVS SSIGTEPDSV EQQFGPGAAR DSQAEPTQPE QAAEAPAEGG PQTNQLETGA SSPERSNEVE MDSTVITVSN SSDKSVLQQP SVNTSIGQPL PSTQLYLRQT ETLTSNLRMP LTLTSNTQVI GTAGNTYLPA LFTTQPAGSG PKPFLFSLPQ PLAGQQQTQFV TVSQPGLSTF TAQLPAPQPL ASSAGHSTAS GQGEKKPYEC TLCNKTFTAK QNYVKHMFVH TGEKPHQCSI CWRSFSLKDY LIKHMVHTTG VRAYQCSICN KRFTQKSSLN VHMRLHRGEK SYECYICKKK FSHKTLLERH VALHSASNGT PPAGTPPGAR AGPPGVVACT EGTTYVCSVC PAKFDQIEQF NDHMRMHVSD G
Source	Mammalian cell
Target Names	ZBTB20
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 20 Alternative name(s): Dendritic-derived BTB/POZ zinc finger protein Zinc finger protein 288
Expression Region	1-741
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.