



Recombinant Human Band 4.1-like protein 4A (EPB41L4A)

Product Code	CSB-MP867204HU
Abbreviation	EPB41L4A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCS5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGCFCAPVEE FYCEVLLLLDE SKLTLTTQQQ GIKKSTKGSV VLDHVFHHVN LVEIDYFGLR YCDRSHQTYW LDPAKTLAEH KELINTGPPY TLYFGIKFYA EDPCKLKEEI TRYQFFLQVK QDVLQGRLPC PVNTAAQLGA YAIQSELGDY DPYKHTAGYV SEYRFVPDQK EELEEAIERI HKTLMGQIPS EAELNYLRTA KSLEMYGVDL HPVYGENKSE YFLGLTPVGV VVYKNKKQVG KYFWPRITKV HFKETQFELR VLGKDCNETS FFFEARSKTA CKHLWKCSVE HHTFFRMPEN ESNSLSRKLS KFGSIRYKHR YSGRTALQMS RDLISQLPRP DQNVTRSRSK TYPKRIAQTQ PAESNSISRI TANMENGNE GTIKIIAPSP VKSFKKAKNE NSPDTQRSKS HAPWEENGPO SGLYNPSDR TKSPKFPYTR RRNPSCGSDN DSVQPVRRRK AHNSGEDSDL KQRRRSR SRC NTSSGSESEN SNREYRKKRN RIRQENDMVD SAPQWEAVLR RQKEKNQADP NNRRSRHRSR SRSPDIQAKE ELWKHIQKEL VDPSGLSEEQ LKEIPYTKIE TQGDPIRIH SHSPRSYRQY RRSQCSDGER SVLSEVNSKT DLVPLPVTR SSDAQSGDA TVHQRRNGSK DSLMEEKPQT STNNLAGKHT AKTIKTIQAS RLKTET
Source	Mammalian cell
Target Names	EPB41L4A
Protein Names	Recommended name: Band 4.1-like protein 4A Alternative name(s): Protein NBL4
Expression Region	1-686
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.