



Recombinant Methanosarcina acetivorans Pyrroline-5-carboxylate reductase (proC)

Product Code	CSB-EP867237MFB
Abbreviation	proC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HH99
Product Type	Recombinant Protein
Immunogen Species	Methanosarcina acetivorans (strain ATCC 35395 / DSM 2834 / JCM 12185 / C2A)
Purity	>85% (SDS-PAGE)
Sequence	MENQKIGFIG AGKMGSALMQ GTIKAGIVTP ENIGASDVYE PFLKDLQAKL GIRVSTDNAV IVRESLIL AVKPQTLSSV LSNLKNEITS EKLVISIAAG VPLSTYEDAL LEGTRVVRVM PNIAATVSEA ASGIAPGKNA TPEDLKAAL IFSAVGTAQ VPESLMDAVT GLSGSGPAFI FPVIEAMADG AVLEGMDRKS ALTAAQTVL GAAKMALETG MHPGELKDMV TSPAGTTIQG IHSLEEAGIR AAFMNAVIRA SERSKELGKK
Source	E.coli
Target Names	proC
Protein Names	Recommended name: Pyrroline-5-carboxylate reductase Short name= P5C reductase Short name= P5CR EC= 1.5.1.2
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.