



# Recombinant Halobacterium salinarum Cell division control protein 6 homolog 6 (cdc6-6)

<b>Product Code</b>	CSB-EP867239HTM-B
<b>Abbreviation</b>	cdc6-6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9HHJ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Halobacterium salinarum (strain ATCC 700922 / JCM 11081 / NRC-1) (Halobacterium halobium)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTPDSSPSSV DDPLFESGHR IFANKDLLKI GHVPEADRIV GRDEEISKLA KRLNGAVHGY SPENVMYIGK TGTGKSLVSK HVCQRAQNAA QDGVEIGTAY IDCAEDNTET QAISLA AKL NDESSTGISV PHTGLSTSKY YKLLWKTLDA QFDSVIIIID EIDLMNDDSV LMKLSRAEEA GKIDCSV GVI AISNKIQYVD NVNERVKSSF QHKELFFKPY DANQLREIMF NREDAFQDGV LSEDVIPLSA AFAAQEHGDA RK AIDILRHA GEVAYEAGAE LVTEEHVRQA QQHA EKDRFR ELVNGAPTQA KAALLATEL SVNSNDDAFL TSRVYDQYEH ICNHLDM DIL SVRRFRDILK EQAFLGVVEI EKINKGSAGG IHLQNRLIED SQVVRETILE DSRMQDWTRE
<b>Source</b>	E.coli
<b>Target Names</b>	orc4
<b>Protein Names</b>	Recommended name: Cell division control protein 6 homolog 6 Short name= CDC6 homolog 6
<b>Expression Region</b>	1-410
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.