



# Recombinant Danio rerio CUGBP Elav-like family member 3 (celf3)

<b>Product Code</b>	CSB-BP867441DIL
<b>Abbreviation</b>	celf3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9IBD1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKEADAIKLF IGQIPRNLEE KDLKPIFEQF GKIYELTVIK DKYTGMHKGC AFLTYCARES ALKAQNALHE QKTLPGMNRP IQVKPADSEG RGDRKLFVGM LGKQLSDADV RKMFEFPGSI EECTVLRGPD GASKGCAFVK YQSNAEAQAA ISALHGSRTL PGASSSLVVK FADTEKERGI RRMQQVASQL GVISPMSLHL GAYNAYTQAL VQQQALVAQS AYLSPVATVA AVQMQQLAAL NPSSIIATPI ASITPSSGTS TPPTIAATPV PALPPPIAVN SYPPVPAAPN GQSASEALYT NGVHPYQAQS PALDPLQQAY TGMQHYTATY PAAYGLVGQP FPHQPTLVAQ QPQQPQQQLQQ REGPEGCNIF IYHLPQEFTD SEMLQMFLPF GNVISAKVVF DRATNQSCKCF GFVSFDNPAS AQAAIQAMNG FQIGMKRLKV QLKRPKDANR PY
<b>Source</b>	Baculovirus
<b>Target Names</b>	celf3
<b>Protein Names</b>	Recommended name: CUGBP Elav-like family member 3 Short name= CELF-3 Alternative name(s): Bruno-like protein 1 CUG-BP- and ETR-3-like factor 3 ELAV-type RNA-binding protein 1 Short name= ETR-1 RNA-binding protein BRUNOL-1
<b>Expression Region</b>	1-452
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.