



Recombinant Mouse RNA binding protein fox-1 homolog 1 (Rbfox1)

Product Code	CSB-MP867494MO
Abbreviation	Rbfox1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JJ43
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNCEREQLRG NQEAAAAPDT MAQPYASAQF APPQNGIPAE YTAPHPHPAP EYTGQTTVPD HTLNLYPPTQ THSEQSADTS AQTVSGTATQ TDDAAPTGGQ PQTQSENTE SKSQPKRLHV SNIPFRFRDP DLRQMFGQFG KILDVEIIFN ERGSKGFGFV TFENSADADR AREKLHGTVV EGRKIEVNNA TARVMTNKKT VNPYTNGWKL NPVVGAVYSP DFYAGTVLLC QANQEGSSMY SGPSSLVYTS AMPGFYPYAA TAAAAYRGAH LRGRGRTVYN TFRAAAPPPP IPAYGGVVYQ DGFYGADIYG GYAAARYAQP TPATAAAYS DSYGRVYAADP YHHTLAPAPT YGVGAMNAFA PLTDAKTRSH ADDVGLVLSS LQASIYRGGY NRFAPY
Source	Mammalian cell
Target Names	Rbfox1
Protein Names	Recommended name: RNA binding protein fox-1 homolog 1 Alternative name(s): Ataxin-2-binding protein 1 Fox-1 homolog A
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.