



Recombinant Mouse Partitioning defective 6 homolog gamma (Pard6g)

Product Code	CSB-EP867510MO-B
Abbreviation	Pard6g
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JK84
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNRSFHKSQT LRFYDCSAVE VKSKFGAEFR RFSLDRHKPG KFEDFYQLVV HTHHISNTEV TIGYADVHGD LLPINDDNF CKAVSSANPL LRVFIQKREE ADHYSFGAGT LSRKKKVLVT LRDDGLRRRP HLNISMPHDF RPVSSIIDVD ILPETHRRVR LYRHGYEKPL GFYIRDGTSV RVTPHGLEKV PGIFISRMVP GGLAESTGLL AVNDEVLEVN GIEVAGKTLQ QVTDMMIANS HNLIVTVKPA NQRNNVVRSS RTSGSSVHST DSTTSHHSLP GAHVLRQNSD VESDEEADIV IEGALEPQHI PKTQAVPPGS LSRANGTSLA HGLHRRDMSL HSSGRESNGS IHRFLSSLKP DPRHSLVLPQ GGVEEHGPAI TL
Source	E.coli
Target Names	Pard6g
Protein Names	Recommended name: Partitioning defective 6 homolog gamma Short name= PAR-6 gamma Short name= PAR6A
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.