



Recombinant Arabidopsis thaliana Axial regulator YABBY 4 (YAB4)

Product Code	CSB-EP867913DOA
Abbreviation	YAB4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LDT3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTKLPNMTT LNHFLDLPQG ICHVQCGFCT TILLVSPFT SLSMVTVRC GHCTSLLSVN LMKASFIPLH LLASLSHLDE TGKEEVAATD GVEEEAWKVN QEKENSPTTL VSSSDNEDED VSRVYQVVK PPEKRQRAPS AYNCFIKKEI RRLKAQNPSM AHKEAFSLAA KNWAHFPPAH NKRAASDQCF CEEDNNAILP CNVFEDHEES NNGFRERKAQ RHSIWGKSPF E
Source	E.coli
Target Names	YAB4
Protein Names	Recommended name: Axial regulator YABBY 4 Alternative name(s): Protein INNER NO OUTER
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.