



Recombinant *Oryza sativa* subsp. japonica Zinc finger protein STAMENLESS 1 (SL1)

Product Code	CSB-MP867932OFG
Abbreviation	SL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LG97
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MNSSRRQEGS PLDLNNLPDE FGKQTVESST TTAASSAEAS RVTKKKSNGG KDEAGKVYEC RFCSLKFCKS QALGGHMNRH RQERETETLN RARQLVFGND SLAAVGAQLN FRDVNMGGGG AAAPPPTMQM GGGGFRGGGV GGDPCIPLRP VQPRLSPPQP PPYHHYLYTT TAPPSALHPM SYPATYPAPP RHQQPAAVGD YVIGHAVSAG DALVAPPPPP HRASFSCFGA PLAAPPANVQ PDNGNCNCSF GCGHSNRNVN AAS
Source	Mammalian cell
Target Names	SL1
Protein Names	Recommended name: Zinc finger protein STAMENLESS 1 Alternative name(s): OsJAG Zinc finger protein OPEN BEAK
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.