



Recombinant *Oryza sativa* subsp. *japonica* Histone H2B.10 (H2B.10)

Product Code	CSB-EP867933OFG-B
Abbreviation	H2B.10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LG12
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	APKAEKKPA AKKPAEEEP AAEKAEKALAG KKPKAERLP AGKAEKSSGE GKKAGRKKAK KSVETYKIYI FKVLKQVHPD IGISSKAMSI MNSFINDIFE KLAGESAKLA RYNKKPTITS REIQTSVRLV LPGELAKHAV SEGTKAVTKF TSA
Source	E.coli
Target Names	H2B.10
Protein Names	Recommended name: Histone H2B.10
Expression Region	2-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.