



Recombinant Arabidopsis thaliana Putative B3 domain-containing protein At3g28852 (At3g28852)

Product Code	CSB-BP867935DOA
Abbreviation	At3g28852
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LH88
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MADNIDLELS LQTNHSMIIN KRLTQSDVDY NNRLHLPKRE FEQFILPEME WELVMNLRNS VEVIVKDVNG NEYHVTLVKY QNGHYFFMGK WMDIVRAKGY KRDDEISLLW DKSNEVFYII
Source	Baculovirus
Target Names	At3g28853
Protein Names	Recommended name: Putative B3 domain-containing protein At3g28852
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.