



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 4 (At1g07360)

Product Code	CSB-BP867979DOA
Abbreviation	At1g07360
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LNV5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAHRILRDHE ADGWERSDFP IICESCLGDN PYVRMTKANY DKECKICTRP FTVFRWRPGR DARYKKTEIC QTCCCKLKNVC QVCLLDLEYG LPVQVRDTAL NISTHDSIPK SDVNREYFAE EHDRKARAGL DYESSFGKMR PNDTILKLQR TTPYYKRNRA HVCSEFFIRGE CTRGAECPPYR HEMPETGELS QQNIKDRYYG VNDPVAMKLL GKAGEMGTLE SPDDESI KTL YVGGLNSRIL EQDIRDQFYA HGEIESIRIL ADKACAFVTY TSREGAEKAA QELSNRLVIN GQRLKLTWGR PKPDQDGANQ QGGVAHSGLL PRAVISQQHN QPPPMQQYYM HPPPANQDKP YYPSMDPQRM GAVISTQEAG GSSTENNGAS SSSYMMPPHQ SYPPPPYGYM PSPYQQQYPP NHHHQPSPMQ HYAPPPAAYP YPQQPGPGSR PAPSPTAVSA ISPDSAPAGS GAPSGSSQQA PDVSTATGSS Q
Source	Baculovirus
Target Names	At1g07360
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 4 Short name= AtC3H4
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.