



Recombinant Arabidopsis thaliana Histone chaperone ASF1B (ASF1B)

Product Code	CSB-MP868003DOA
Abbreviation	ASF1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LS09
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSSINITNVT VLDNPAPFVN PFQFEISYEC LTSLKDDLEW KLIYVGS AED ETYDQVLESV LVGPVNVGN Y RFVLQADSPD PLKIREEDII GVTVLLLTCS YMDQEFIRVG YYVNNDYDDE QLREEPPTKV LIDKVQRN IL TDKPRVTKFP INFHPENEQT LGDGPAPTEP FADSVVNGEA PVFLEQPQKL QEIEQFDDSD VNGEIALLD QPQNLQET
Source	Mammalian cell
Target Names	ASF1B
Protein Names	Recommended name: Histone chaperone ASF1B Alternative name(s): Anti-silencing function protein 1-like protein b Short name= Anti-silencing function 1b protein Silencing group A protein 1
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.