



# Recombinant Arabidopsis thaliana Geranylgeranyl pyrophosphate synthase 9, chloroplastic (GGPPS9)

<b>Product Code</b>	CSB-EP868030DOA
<b>Abbreviation</b>	GGPPS9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LUE1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	S ALTSQDAGHM IQPEGKSNNDN NSAFDFKLYM IRKAESVNAA LDVSVPLLKP LTIQEAVRYS LLAGGKRVRP LLCIAACELV GGDEATAMSA ACAVEMIHTS SLIHDDLPCM DNADLRRGKP TNHVKYGEDM AVLAGDALLA LAFEHMTVVS SGLVAPEKMI RAVVELARAI GTTGLVAGQM IDLASERLNP DKVGGLEHLEF IHLHKTAALL EAAAVLGVIM GGGTEQEIEK LRKYARCIGL LFQVVDDILD VTKSTEELGK TAGKDV MAGK LTYPRLIGLE GSREVAEKLR REAEEQLLGF DPSKAAPLVA LASYIACRHN
<b>Source</b>	E.coli
<b>Target Names</b>	GGPPS9
<b>Protein Names</b>	Recommended name: Geranylgeranyl pyrophosphate synthase 9, chloroplastic Short name= GGPP synthase 9 Short name= GGPS9 EC= 2.5.1.- Alternative name(s): (2E,6E)-farnesyl diphosphate synthase 9 Dimethylallyltranstransferase
<b>Expression Region</b>	40-360
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.