



Recombinant Arabidopsis thaliana Putative pirin-like protein At3g59260 (At3g59260)

Product Code	CSB-EP868058DOA-B
Abbreviation	At3g59260
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LX45
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKVMVLSLGK ASPSKSDHEL LDPFVSLVEF SVSPPPGGFKD HPHRGFESVT YMFQGGIIHQ DCNGNKGTH EGDVQWMTAG RGIHSEMPE EQVNKGLQLW INLPSSAKMI EPKNIEISSS EIPSADDYGV EVKVIAGESM GVKSPFYTKT PIMFLDFTLD PKAQTHQAVP ESWTAFAYIV EGDEGVFSSS DSSTVQAHNV VVFVGTGDEVS VVWNTSNRPL RFLIAGEPI GEPVVQHGPF VMNSQDEIEM TIGDYRNGMN GFEMAKHWRS E
Source	E.coli
Target Names	At3g59260
Protein Names	Recommended name: Putative pirin-like protein At3g59260
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.