



# Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At5g02910 (At5g02910)

<b>Product Code</b>	CSB-MP868078DOA
<b>Abbreviation</b>	At5g02910
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LYZ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MITGDSIDGM DFISSLPDEI LHHILSSVPT KSAIRTSLLS KRWRYVWSET PSLSIDCRRRA DPNSIDKTLS FFSAPKITSF HLHTLLNRI DSVNGCIEFA ISHNAEKLSL ESRDYRVRNY KFPDFFYTNS SVKQLFVDSG SVHLIPRCTV SWTSLKNLSL SNCTLSDES F LKILSGSPLL ESLELLYCAE YMCLDLSQSQ HLRRLEIDRS DWFMGPTKIV APHLHCLRLR HSRLPCSLVD VSSLTEADLN IYFGDLKTVT AGFLQHN VVK MLQMLQNVEK LTIGGTFLQM LSLAALCGVP FPTLKVKT LT LETMIIRSVI PGITKLLRYT PGLRKLTIHT VKCSSISELH LNDYLRKHSL NQRQCWRSKD SVFPGSLETI SMLASKHAES NLVALFMERL LKSTKSLETM VVLLVDYLDA LGFEELLAMA TTLSHNNDVS VLIKRSYIKY VSNTFPQR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At5g02910
<b>Protein Names</b>	Recommended name: F-box/LRR-repeat protein At5g02910
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.