



Recombinant Arabidopsis thaliana Homocysteine S-methyltransferase 2 (HMT-2)

Product Code	CSB-EP868108DOA-B
Abbreviation	HMT-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M1W4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTGNSFNSMK DFLKQTGGYA VIDGGLATEF ERHGADLNDP LWSAKCLVTS PHLIHTVHLD YLEAGADIIS SASYQATIQQ FEAKEGFSREE SESLLKKSVE IATEARNSYY DKCGTSSSMD DKILKKRPIL VAASVGSYGA YLADGSEYSG IYGDSITLEK LKDFHRRRLQ VLAESGADLI AFETIPNKIE AQAFADLLEE GDVKIPGWFS FNSKDG VNVV SGDSIKECIS IAENCEKVVA VGINCTPPRF IEGLVLEIEK VTSKIPILVYP NSGESYDADR KEWVENTGVG DEDFVSYVEK WMDAGVSLLG GCCRTTPTTI RAIHKRLVNR RSL
Source	E.coli
Target Names	HMT-2
Protein Names	Recommended name: Homocysteine S-methyltransferase 2 EC= 2.1.1.10 Alternative name(s): S-methylmethionine:homocysteine methyltransferase 2 Short name= AtHMT-2 Short name= SMM:Hcy S-methyltransferase 2
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.