



Recombinant Arabidopsis thaliana Probable mannose-1-phosphate guanylyltransferase 2 (At3g55590)

| | |
|--------------------------|--|
| Product Code | CSB-EP868115DOA |
| Abbreviation | At3g55590 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9M2S0 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MKALILVGGF GTRLRPLTLS LPKPLVDFAN KPMILHQIEA LKAIGVDEVV LAINYEPEQL LVMSKFSNDV EATLGIKITC SQETEPLGTA GPLALARDKL VDGSGQPFFV LNSDVISDYP LEEMIAFHNA HGGEASIMVT KVDEPSKYGV VVMEEATGRV ERFVEKPKLF VGNKINAGIY LLNPSVLDR I ELRPTSIEKE IFPQIAEAEK LYAMLLPGFW MDIGQPRDYI TGLRLYLDSL RKKSPSKLAT GPHILGNVLV DETAEIGEGC LIGPNVAIGP GCVVESGVRL SHCTVMRGVH VKRYACISSS IIGWHSTVGQ WARVENMSIL GKNVYVCDEI YCNGGVVLHN KEIKSDILKP DIVM |
| Source | E.coli |
| Target Names | At3g55590 |
| Protein Names | Recommended name: Probable mannose-1-phosphate guanylyltransferase 2 EC= 2.7.7.13 |
| Expression Region | 1-364 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.