



Recombinant Bovine Plasma serine protease inhibitor (SERPINA5)

Product Code	CSB-EP868236BO
Abbreviation	SERPINA5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9N2I2
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	KKIQEV PPAVTTAPPG SRDFVFDLYR ALAAAAPAQN IFFSPLSITV SLAMLSLGAQ SNTKAQILEG LGIGPGEGSE EELHSASQRL LRELQQPQDS LQLSLGNALF TKPRLPIQEA FLGAMRTLIL ADTFPTNFED PEGAKKKIND YVAKQTKGKI VDLIKSLDGT QVMVMVNYIF FKAKWETSFN LKSTHEQDFY VTPETVVRVP MMKQQDQFYY LLDRNLCKV VGVVYQGNAT AFFILPREGE MEQVENGLKE KTLKKWLRMP MKRRLELYLP KFSIEGSYQL EEVLPKLGIR DIFTSDADLT GISNHSSIRV SEMVHKAVVE VDESGTQAAA ATGMVITFKS ARLGSQRIVF NRPFLVLIVK NSKHILFLGK VTRP
Source	E.coli
Target Names	SERPINA5
Protein Names	Recommended name: Plasma serine protease inhibitor Alternative name(s): Protein C inhibitor Short name= PCI Serpin A5
Expression Region	25-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.