



Recombinant Human Lysine-rich coiled-coil protein 1 (KRCC1)

Product Code	CSB-EP868279HU
Abbreviation	KRCC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPI7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKHSKPTYDS FQDELEDYIK VQKARGLEPK TCFRKMKGDY LETCGYKGEV NSRPTYRMFD QRLPSETIQT YPRSCNIPQT VENRLPQWLP AHDSRLRLDS LSYCQFTRDC FSEKPVPLNF NQQEYICGSH GVEHRVYKHF SSDNSTSTHQ ASHKQIHQKR KRHPPEGREK SEEERSKHKR KKSCEEIDLD KHKSIQRKKT EVEIETVHVS TEKLKNRKEK KSRDVVSKKE ERKRTKKKKE QGQERTEEEM LWDQSILGF
Source	E.coli
Target Names	KRCC1
Protein Names	Recommended name: Lysine-rich coiled-coil protein 1 Alternative name(s): Cryptogenic hepatitis-binding protein 2
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.