



Recombinant Human Transcriptional repressor scratch 2 (SCRT2)

Product Code	CSB-BP868281HU
Abbreviation	SCRT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQ03
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRSFLVKKI KGDGFQCSGV PPTYHPLET AYVLPGARGP PGDNGYAPHR LPPSSYDADQ KPGLELAPAE PAYPPAAPEE YSDPESPQSS LSARYFRGEA AVTDSYSMDA FFISDGRSRR RRGGGGGDAG GSGDAGGAGG RAGRAGAQAG GGHRHACAEC GKTYATSSNL SRHKQTHRSL DSQLARKCPT CGKAYVSMIPA LAMHLLTHNL RHKCGVCGKA FSRPWLLQGH MRSHTGEKPF GCAHCGKAFA DRSNLRAHMQ THSAFKHYRC RQCDKSFALK SYLHKHCEAA CAKAAEPPPP TPAGPAS
Source	Baculovirus
Target Names	SCRT2
Protein Names	Recommended name: Transcriptional repressor scratch 2 Alternative name(s): Scratch homolog 2 zinc finger protein
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.