



Recombinant Human Nuclear receptor-interacting protein 3 (NRIP3)

Product Code	CSB-MP868283HU
Abbreviation	NRIP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQ35
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFYSGLLTEG GRKETDMREA ASLRQQRRMK QAVQFIHKDS ADLLPLDGLK KLGSSKDMQP HNILQRRRLME TNLSKLRSGP RVPWASKTNK LNQAKSEGLK KSEEDDMILV SCQCAGKDVK ALVDTGCLYN LISLACVDRL GLKEHVKSHK HEGEKLSLPR HLKVVGQIEH LVITLGLSRL DCPAAVVDDN EKNLSLGLQT LRSLKCIINL DKHRLIMGKT DKEEIPFVET VSLNEDNTSE A
Source	Mammalian cell
Target Names	NRIP3
Protein Names	Recommended name: Nuclear receptor-interacting protein 3 Alternative name(s): Sarcoma antigen NY-SAR-105
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.