



Recombinant Human SNF-related serine/threonine-protein kinase (SNRK)

| | |
|--------------------------|---|
| Product Code | CSB-EP868301HU-B |
| Abbreviation | SNRK |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9NRH2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAGFKRGYDG KIAGLYDLDK TLGRGHFAVV KLARHVFTGE KVAVKVIDKT KLDTLATGHL FQEVRCMKLV QHPNIVRLYE VIDTQTKLYL ILELGDGGDM FDYIMKHEEG LNEIDLAKKYF AQIVHAISYC HKLHVVHRDL KPENVVFFEK QGLVKLTDFG FSNKFQPGKK LTTSCGSLAY SAPEILLGDE YDAPAVDIWS LGVILFMLVC GQPPFQEAND SETLTMIMDC KYTVPSHVSK ECKDLITRML QRDPKRRASL EEIENHPWLQ GVDPSPATKY NIPLVSYKNL SEEEHNSIIQ RMVLGDIADR DAIVEALETN RYNHITATYF LLAERILREK QEKEIQTRSA SPSNIKAQFR QSWPTKIDVP QDLEDDLTAT PLSHATVPQS PARAADSVLN GHRSKGLCDS AKKDDLPELA GPALSTVPPA SLKPTASGRK CLFRVEEDEE EDEEDKKPMS LSTQVVLRRK PSVTNRLTSR KSAPVLNQIF EEGESDDEFD MDENLPPKLS RLKMNIASPG TVHKRYHRRK SQGRGSSCSS SETSDDDDSES RRRLDKDSGF TYSWHRRDSS EGPPGSEGDG GGQSKPSNAS GGVDKASPSE NNAGGGSPSS GSGGNPTNTS GTTRRCAGPS NSMQLASRSA GELVESLKLM SLCLGSQLHG STKYIIDPQN GLSFSSVKVQ EKSTWKM CIS STGNAGQVPA VGGIKFFSDH MADTTTELER IKSKNLKNNV LQLPLCEKTI SVNIQRNPKE GLLCASSPAS CCHVI |
| Source | E.coli |
| Target Names | SNRK |
| Protein Names | Recommended name: SNF-related serine/threonine-protein kinase EC=2.7.11.1 Alternative name(s): SNF1-related kinase |
| Expression Region | 1-765 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.