



Recombinant Human SNF-related serine/threonine-protein kinase (SNRK)

Product Code	CSB-EP868301HU
Abbreviation	SNRK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRH2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGFKRGYDG KIAGLYDLDK TLGRGHFAVV KLARHVFTGE KVAVKVIDKT KLDTLATGHL FQEVRCMKLV QHPNIVRLYE VIDTQTKLYL ILELGDGGDM FDYIMKHEEG LNEIDLAKKYF AQIVHAISYC HKLHVVHRDL KPENVVFFEK QGLVKLTDFG FSNKFQPGKK LTTSCGSLAY SAPEILLGDE YDAPAVDIWS LGVILFMLVC GQPPFQEAND SETLTMIMDC KYTVPSHVSK ECKDLITRML QRDPKRRASL EEIENHPWLQ GVDPSPATKY NIPLVSYKNL SEEEHNSIIQ RMVLGDIADR DAIVEALETN RYNHITATYF LLAERILREK QEKEIQTRSA SPSNIKAQFR QSWPTKIDVP QDLEDDLTAT PLSHATVPQS PARAADSVLN GHRSKGLCDS AKKDDLPELA GPALSTVPPA SLKPTASGRK CLFRVEEDEE EDEEDKKPMS LSTQVVLRRK PSVTNRLTSR KSAPVLNQIF EEGESDDEFD MDENLPPKLS RLKMNIASPG TVHKRYHRRK SQGRGSSCSS SETSDDDDSES RRRLDKDSGF TYSWHRRDSS EGPPGSEGDG GGQSKPSNAS GGVDKASPSE NNAGGGSPSS GSGGNPTNTS GTTRRCAGPS NSMQLASRSA GELVESLKLM SLCLGSQLHG STKYIIDPQN GLSFSSVKVQ EKSTWKMCI S STGNAGQVPA VGGIKFFSDH MADTTTELER IKSKNLKNNV LQLPLCEKTI SVNIQRNPKE GLLCASSPAS CCHVI
Source	E.coli
Target Names	SNRK
Protein Names	Recommended name: SNF-related serine/threonine-protein kinase EC= 2.7.11.1 Alternative name(s): SNF1-related kinase
Expression Region	1-765
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.