



Recombinant Human SPATS2-like protein (SPATS2L)

Product Code	CSB-MP868324HU
Abbreviation	SPATS2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NUQ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AELNTHVNV KEKIYAVRSV VPNKSNN EIV LVLQQFDFNV DKAVQAFVDG SAIQVLKEWN MTGKKKNNKR KRSKSKQHQQ NKDAKDKVER PEAGPLQPQP PQIQNGPMNG CEKDSSTDS ANEKPALIPR EKKISILEEP SKALRGVTEG NRLLQQKLSL DGNPKPIHGT TERSDGLQWS AEQPCNPSKP KAKTSPVKS N TPAAHLEIKP DELAKKRGPN IEKSVKDLQR CTVSLTRYRV MIKEEVDSSV KKIKAAFAEL HNCIIDKEVS LMAEMDKVKE EAMEILTARQ KKAEEKRLT DLASQMAEMQ LAELRAEIKH FVSEKRYDEE LGKAARFSCD IEQLKAQIML CGEITHPKNN YSSRTPCSSL LPLLNAHAAT SGKQSNFSRK SSTHNKPSEG KAAPKMOVSS LPSTADPSHQ TMPANKQNGS SNQRRRFNPQ YHNNRLNGPA KSQGS GNEAE PLGKGNSRHE HRRQPHNGFR PKNKGGAKNQ EASLGMKTPE APAHSEKPRR RQHAADTSEA RPFGRSVGRV SQCNLCPTRI EVSTDAAVLS VPAVTLVA
Source	Mammalian cell
Target Names	SPATS2L
Protein Names	Recommended name: SPATS2-like protein Alternative name(s): DNA polymerase-transactivated protein 6 Stress granule and nucleolar protein Short name= SGNP
Expression Region	2-558
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.