



# Recombinant Human SPATS2-like protein (SPATS2L)

<b>Product Code</b>	CSB-YP868324HU
<b>Abbreviation</b>	SPATS2L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NUQ6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AELNTHVNV KEKIYAVRSV VPNKSNN EIV LVLQQDFDNV DKAVQAFVDG SAIQVLKEWN MTGKKKNNKR KRSKSKQHGG NKDAKDKVER PEAGPLQPQP PQIQNGPMNG CEKDSSTDS ANEKPALIPR EKKISILEEP SKALRGVTEG NRLLQQKLSL DGNPKPIHGT TERSDGLQWS AEQPCNPSKP KAKTSPVKS NTPAAHLEIKP DELAKKRGPN IEKSVKDLQR CTVSLTRYRV MIKEEVDSSV KKIKAAFAEL HNCIIDKEVS LMAEMDKVKE EAMEILTARQ KKAEEKRLT DLASQMAEMQ LAELRAEIKH FVSEKRYDEE LGKAARFSCD IEQLKAQIML CGEITHPKNN YSSRTPCSSL LPLLNAHAAT SGKQSNFSRK SSTHNKPSEG KAAPKMOVSS LPSTADPSHQ TMPANKQNGS SNQRRRFNPQ YHNNRLNGPA KSQGS GNEAE PLGKGNSRHE HRRQPHNGFR PKNKGGAKNQ EASLGMKTPE APAHSEKPRR RQHAADTSEA RPFGRSVGRV SQCNLCPTRI EVSTDAAVLS VPAVTLVA
<b>Source</b>	Yeast
<b>Target Names</b>	SPATS2L
<b>Protein Names</b>	Recommended name: SPATS2-like protein Alternative name(s): DNA polymerase-transactivated protein 6 Stress granule and nucleolar protein Short name= SGNP
<b>Expression Region</b>	2-558
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .