



# Recombinant Human Protein asunder homolog (Asun)

<b>Product Code</b>	CSB-EP868331HU
<b>Abbreviation</b>	Asun
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NVM9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKIFSESHKT VFVVDHCPYM AESCRQHVEF DMLVKNRTQG IIP LAPISKS LWTCVSESSM EYCRIMYDIF PFKKLVNFIV SDSGAHLVNS WTQEDQNLQE LMAALAAVGP PNPRADPECC SILHGLVAAV ETLCKITEYQ HEARTLLMEN AERVGNRGRI ICITNAKSDS HVRMLEDCVQ ETIHEHNKLA ANSDHLMQIQ KCELVLIHTY PVGEDSLVSD RSKKELSPVL TSEVHSVVRAG RHLATKLNIL VQQHFDLAST TITNIPMKEE QHANTSANYD VELLHHKDAH VDFLKSGDSH LGGGSREGSF KETITLKWCT PRTNNIELHY CTGAYRISPV DVNSRPSSCL TNFLLNGRSV LLEQPRKSGS KVISHMLSSH GGEIFLHVLS SRSILEDPP SISEGCGGRV TDYRITDFGE FMRENRLTPF LDPYKIDGS LEVPLERAKD QLEKHTRYWP MIISQTTIFN MQAVVPLASV IVKESLTEED VLNCQKTIYN LVDMERKNDP LPISTVGTRG KGPKRDEQYR IMWNELETLV RAHINNSEKH QRVLECLMAC RSKPPEEEER KKRGRKREDK EDKSEKAVKD YEQEKSWQDS ERLKGILERG KEELAEAEII KDSPDSPEPP NKKPLVEMDE TPQVEKSKGP VSLLSLWSNR INTANSRKHQ EFAGRLNSVN NRAELYQHLLK EENGMETTEN GKASRQ
<b>Source</b>	E.coli
<b>Target Names</b>	INTS13
<b>Protein Names</b>	Recommended name: Protein asunder homolog Alternative name(s): Cell cycle regulator Mat89Bb homolog Sarcoma antigen NY-SAR-95
<b>Expression Region</b>	1-706
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.