



Recombinant Human Protein asunder homolog (Asun)

Product Code	CSB-MP868331HU
Abbreviation	Asun
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NVM9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKIFSESHKT VFVVDHCPYM AESCRQHVEF DMLVKNRTQG IIP LAPISKS LWTCSVESSM EYCRIMYDIF PFKKLVNFIV SDSGAHLVNS WTQEDQNLQE LMAALAAVGP PNPRADPECC SILHGLVAAV ETLCKITEYQ HEARTLLMEN AERVGNRGR IICITNAKSDS HVRMLEDCVQ ETIHEHNKLA ANSDHLMQIQ KCELVLIHTY PVGEDSLVSD RSKKELSPVL TSEVHSVVRAG RHLATKLNIL VQQHFDLAST TITNIPMKEE QHANTSANYD VELLHHKDAH VDFLKSGDSH LGGGSREGSF KETITLKWCT PRTNNIELHY CTGAYRISPV DVNSRPSSCL TNFLLNGRSV LLEQPRKSGS KVISHMLSSH GGEIFLHVLS SRSILEDPP SISEGCGGRV TDYRITDFGE FMRENRLTPF LDPYKIDGS LEVPLERAKD QLEKHTRYWP MIISQTTIFN MQAVVPLASV IVKESLTEED VLNCQKTIYN LVDMERKNDP LPISTVGTRG KGPKRDEQYR IMWNELETLV RAHINNSEKH QRVLECLMAC RSKPPEEEER KKRGRKREDK EDKSEKAVKD YEQEKSWQDS ERLKGILERG KEELAEAEII KDSPDSPEPP NKKPLVEMDE TPQVEKSKGP VSLLSLWSNR INTANSRKHQ EFAGRLNSVN NRAELYQHLLK EENGMETTEN GKASRQ
Source	Mammalian cell
Target Names	INTS13
Protein Names	Recommended name: Protein asunder homolog Alternative name(s): Cell cycle regulator Mat89Bb homolog Sarcoma antigen NY-SAR-95
Expression Region	1-706
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.