



Recombinant Human PCNA-interacting partner (PARPBP)

Product Code	CSB-EP868339HU-B
Abbreviation	PARPBP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWS1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAVFNQKSVS DMIKEFRKNW RALCNSERTT LCGADSMMLLA LQLSMAENNK QHSGEFTVSL SDVLLTWKYL LHEKLNLPVE NMDVTDHYED VRKIYDDFLK NSNMLDLIDV YQKCRALTSN CENYNTVSPS QLLDFLSGKQ YAVGDETDLS IPTSPTSKYN RDNEKVQLLA RKIIFSYNL LVNSKNDLAV AYILNIPDRG LGREAFTDLK HAAREKQMSI FLVATSFIRT IELGGKGYAP PPSDPLRTHV KGLSNFINFI DKLDEILGEI PNPSIAGGQI LSVIKMQLIK GQNSRDPFCK AIEEVAQDLD LRIKNIINSQ EGVVALSTTD ISPARPKSHA INHGTAFCGR DTVKALLVLL DEEAANAPTK NKAELLYDEE NTIHHHGTSI LTLFRSPTQV NNSIKPLRER ICVSMQEKKI KMKQTLIRSQ FACTYKDDYM ISKDNWNNVN LASKPLCVLY MENDLSEGVN PSVGRSTIGT SFGNVHLDRS KNEKVSRSKST SQTGNKSSKR KQVDLDGENI LCDNRNEPPQ HKNAKIPKKS NDSQNRLYGK LAKVAKSNKC TAKDKLISGQ AKLTQFFRL
Source	E.coli
Target Names	PARPBP
Protein Names	Recommended name: PCNA-interacting partner Short name= PARI Alternative name(s): PARP-1 binding protein PARP1-binding protein Short name= PARPBP
Expression Region	1-579
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.